


# MedArt FRx Intenz CO<sub>2</sub>

## MedArt Solutions for Medical & Aesthetic Indications:

- FRx Skin Rejuvenation
- Wrinkles
- Warts
- Age Spots
- Scars
- Stretch Marks
- Tissue Ablation
- Laser Vaginal Rejuvenation  
and many more ..



MedArt FRx Intenz CO<sub>2</sub> 10600 nm is a highly versatile and effective system for a full ablation and surgical treatments. In addition, the advanced unique MedArt 3D scanner enables a full fractional treatment. MedArt fractional treatment can be simply individualized by selecting density, pattern and size of treatment area. The compact, portable, ergonomic design of MedArt Intenz with articulated arm & scanner makes it comfortable and easy to use.



# MedArt

# MedArt

## MedArt FRx Intenz CO<sub>2</sub> Specification

Laser Type:	Continuous Wave Carbon Dioxide (CO <sub>2</sub> ) Laser, Class 4 (EEC 60825)
Output Power Range:	0.1-25 W in steps of 0.1W
Wavelength:	10.600 nm
Aiming Beam:	Adjustable Laser Diode (635 nm), 2 x Crossing Beams
Mode of Operation:	Continuous or Pulsed Operation
Frequency Range:	0.3 - 100 Hz
Pulse Width:	10 – 990 msec
Beam Delivery System:	Articulated Arm , 7 mirrors
Arm length:	165 cm
Hand Piece:	Fixed Focus (0.165 mm), Adjustable Variospot (0.165 – 1.00 mm), Scanner
Spot Size:	0.165 – 1 mm
Main Voltage:	100 – 240 V, 50/60 Hz
Safety Class:	Electrical Class I Type BF
Size:	58 x 46 x 17 cm
Weight:	Approx. 21 kg including Articulated Arm
Product Code:	MedArt 620
Certifications:	 0459 / FDA
<b>Scanner Specifications:</b>	
Scan Size:	15x15 or 6x6 mm
Scan Density:	Low 3% - Medium 7% - Max 12%
Scan Dwell Time:	0.5, 1, 1.5, 2, 5, 7 msec
Treatment Distance:	48 mm
Product Code:	MedArt 459
Certifications:	 0459 / FDA

MedArt A/S is a Danish company with an international profile developing, manufacturing, and selling world class medical laser for aesthetic and medical dermatology. The company was founded in 1979 and the head office, R&D, manufacturing, and marketing department are all located in Copenhagen, Denmark. We offer a broad range of medical lasers combined with the MedArt range of scanners and MedArt cooling systems for most dermatological indications.



**Denmark**  
**Industriholmen 15A**  
**2650 Hvidovre**  
**Phone: +45 36 34 23 00**  
**E-mail: [info@medart.dk](mailto:info@medart.dk)**  
**Website: [www.medart.dk](http://www.medart.dk)**

**Germany**  
**Rudendorfer Weg 39**  
**D 96188 Stettfeld**  
**Phone: +49 (0)9522 707222**

**USA**  
**1850 Diamond Street No 104**  
**San Marcos, CA 92078**  
**Phone: 1-760-798-9550**



**MedArt VariMed™** The unique compact 810nm Diode Laser System in combination with 3D-scanner and cooling system. For treatments of various vascular lesions, transcutaneous leg vein treatment, endovenous veins treatment as well as hair removal on all body areas.



**MedArt SmartSculpt™** The compact, table-top sized, safe, and accurate 980nm diode laser system for assisted lipolysis, hyperhidrosis, body contouring and sculpting.



**MedArt FRx CO<sub>2</sub> Laser** The well-proven CO<sub>2</sub> platform to offer advanced fractional treatments with the unique MedArt scanner which distributing the energy with the adjusted pulse width and spot density. The scanner is connected to the laser by a unique flexible optical fiber.



**MedArt ILVO™ 1470nm Diode Laser** The compact and advanced 1470nm Diode Laser system to maximize the comfort, safe, fast and effectiveness of the varicose vein treatments. MedArt ILVO can also be supplied in 980nm wavelength in request of the clinics.